

Tactical Hemorrhage Control Trainer (THCT)

Operative | Experience

Copyright Notice

Copyright ©2018 Operative Experience, Inc. All rights reserved.

This documentation is protected by the copyright laws as an unpublished work.

Confidentiality and Proprietary Rights

This documentation is the confidential property of Operative Experience, Inc. It is furnished under an agreement with Operative Experience, Inc. and may only be used in accordance with the terms of that agreement. The use of this documentation is restricted to customers of Operative Experience, Inc. and their employees who have a need to know the contents of this documentation. The user of this documentation agrees to protect the Operative Experience, Inc. proprietary rights as expressed herein. The user further agrees not to permit access to this documentation by any person for any purpose other than as an aid in the use of the associated products. In no case will this documentation be examined for the purpose of copying any portion of the product described herein or to design another product to accomplish similar results. This documentation or portions of it may not be copied, modified or distributed without written permission from Operative Experience, Inc. The information in this documentation is subject to change without notice.

No part of this documentation may be reproduced in any form or incorporated into any information retrieval system, electronic or mechanical, without the written permission of the copyright owner. Inquiries regarding permission for use of material contained in this documentation should be addressed to:

Operative Experience, Inc.
500 Principio Parkway West, Suite 900
North East, MD 21901
Phone: 410-287-5748
E-mail: info@operativeexperience.com

Disclaimer

This documentation is provided “as is” without any express or implied warranty of any kind including warranties of merchantability, non-infringement of intellectual property, or fitness for any particular purpose. In no event shall Operative Experience, Inc. or its representatives be liable for any damages whatsoever (including, without limitation, damages for loss of profits, business interruption, or loss of information) arising out of the use of or inability to use this documentation, even if Operative Experience, Inc. has been advised of the possibility of such damages. Because some jurisdictions prohibit the exclusion or limitation of liability for consequential or incidental damages, the above limitation may not apply to you. Operative Experience, Inc. further does not warrant the accuracy or completeness of the information, text, graphics, or other items contained within these materials.

Operative Experience is a registered trademark of Operative Experience, Inc. Soldier Strong is a trademark of Operative Experience, Inc. Velcro is a registered trademark of Velcro BVBA. Germaben II is a registered trademark of International Specialty Products. All other company names, logos, and trademarks mentioned herein are the property of their respective owners.

Table of Contents

Figures.....	v
Tactical Hemorrhage Control Trainer	1
Product Components	1
External Access	1
Caring for Your Simulator.....	2
Lubrication	2
Use and Treatment	2
Cleaning.....	2
Storage	2
Troubleshooting.....	2
Understanding Your Simulator	3
Common Features.....	3
Operating Your Simulator	4
Adding Blood to Your Simulator	4
Using the Key Fob to Control the Bleeding.....	5
Shutting Down the Simulator.....	5
Removing the Blood Mixture	5
Flushing the System	6
Storing after Use.....	6
Performing Interventions.....	6
Needle Decompression.....	6
Oropharyngeal Airway	7
Nasopharyngeal Airway	7
Wound Packing	8
Tourniquet Application	8
Maintaining Your Simulator	9
Emptying/Cleaning the Bleeding System.....	9
Repairing the Skin	9
Charging/Replacing the Battery Pack	10

Removing and Installing the Battery Pack 10

Charging the Battery Pack..... 10

Storing the Simulator 11

Storage and Transportation Instructions..... 11

Using the Storage Base for Transport 11

Figures

FIGURE 1: THE TACTICAL HEMORRHAGE CONTROL TRAINER	3
FIGURE 2: BLOOD INFUSION PORT	4
FIGURE 3: THE HAND PUMP	4
FIGURE 4: USING THE KEY FOB.....	5
FIGURE 5: NEEDLE DECOMPRESSION	6
FIGURE 6: OROPHARYNGEAL AIRWAY	7
FIGURE 7: NASOPHARYNGEAL AIRWAY	7
FIGURE 8: WOUND PACKING	8
FIGURE 9: DIMPLED SKIN FROM NEEDLE PUNCTURES	9
FIGURE 10: STRETCHING THE PUNCTURED SKIN.....	9
FIGURE 11: APPLYING FILLER.....	9
FIGURE 12: RUBBING FILLER INTO THE HOLES	9
FIGURE 13: GURNEY-LIKE STORAGE BASE USED FOR TRANSPORT	11

Tactical Hemorrhage Control Trainer

The Tactical Hemorrhage Control Trainer is a full-sized, life-accurate, remotely-activated simulator designed for point-of-injury, tactical medicine. Ultra-realistic and fully mobile, it is ideally suited for training rapid assessment and treatment of trauma injuries associated with active shooter, disaster and terrorist incidents.

Product Components

Each simulator comes with all the components needed to perform the associated trauma scenarios.

Case with Casters	1
Humeral and sternal I/O replacements	10
Battery Pack NIMH 12V 3.8Ah 45.6 (Wh)	2
Battery Charger	1
Concentrated Fake Blood (1 Gallon)	1
Hand Pump/Canister	1
Leg Bands	3
Neck bands	3
Repair Kits	2
Key Fob Bleeding Controller	1
Gym Shorts	1

External Access

The simulators are designed for realism, with only a few seams or external access points:

- A battery pack access compartment in the upper, inner left thigh
- A blood infusion port located on the underside of the penis

Caring for Your Simulator

Soldier Strong™ but not indestructible! Treat the simulator as you would a combat casualty.

Technical Support

Simulators must be serviced by qualified OEI technicians only.

To schedule service, contact:

Operative Experience Inc.
500 Principio Parkway West
Suite 900
North East, MD 21901

Phone: 410-287-5748

Email: support@operativeexperience.com

Lubrication

Glycerin or mineral oil should be used liberally for enabling the oropharyngeal, nasopharyngeal, and the emergency (surgical) airway on all models.

Use and Treatment

The simulators are durable and Soldier Strong™, but not indestructible. Boots or heel cups and pants/shorts must be worn to protect the skin if the body is dragged during field use. Tissues and joints are designed to replicate human patients. Proper use and treatment promotes proper long-term functionality. The simulators are water-resistant, but not waterproof.

Cleaning

Simulators can be wiped clean with a simple solution of mild soap and water.

Storage

Simulators should be dry, fluid free, and without clothes during long-term storage. See *Storing the Simulator* below. The blood flow system does not need to be flushed for overnight storage but must be for long-term (one week or more) storage. See *Emptying/Cleaning the Bleeding System* below.

Troubleshooting

User-troubleshooting is restricted to the functions described in this manual. Simulator electronics are inaccessible and should only be serviced by qualified OEI technicians. A bleeding failure is generally resolved by charging or replacing the battery. See *Charging/Replacing the Battery Pack* below. For minor tears in the skin, see *Repairing the Skin* below.

Understanding Your Simulator

The Tactical Hemorrhage Control Trainer is a full-sized, ruggedized, soldier-form remotely-activated simulator for rapid assessment and point-of-injury hemorrhage control featuring multiple injuries, including gunshot wounds, stab wounds and crushing injuries.



Figure 1: The Tactical Hemorrhage Control Trainer

Common Features

The Tactical Hemorrhage Control Trainer has the following features:

- Wireless Remote-Controlled
- Needle Decompression
- Oropharyngeal Airway
- Nasopharyngeal Airway
- Wound Packing
- Tourniquet Application

Operating Your Simulator

Prepare your simulator by filling the reservoir of the hand pump with blood, and then use the Controller to prime and operate the simulator.

Adding Blood to Your Simulator

Before operation, infuse the simulator with blood using the hand pump. A gallon of concentrated fake blood is included with the simulator with instructions on its use. Contact Operative Experience Customer Support (support@operativeexperience.com) for additional blood supplies.

To infuse the simulator with blood:

1. Turn the simulator off.

If the simulator has been in use, stop all blood flow prior to infusing, or pressure lock may occur and blood hose coupling will not engage.

2. Fill the reservoir on the hand pump with simulated blood.
3. Lift the penis and extend the blood fill port.
4. Attach the blood fill hose by snapping onto the blood fill port on the simulator located under the penis. (Be sure the hose nozzle is free of debris.)
5. Retract the handle and pump periodically until back-flow check valve on handle emits fluid. This will be an indicator that your reservoir is full.



Figure 2: Blood Infusion Port



Figure 3: The Hand Pump

Using the Key Fob to Control the Bleeding

Your simulator comes equipped with a remote control “key fob” controller to turn on/off the non-pulsatile bleeding system.


To control the bleeding with the controller:

Power on the simulator by pressing the latching illuminated switch on the back of the neck. A blue light will appear on the back of the simulator's neck when powered on.

2. Infuse the blood using the steps described in above.

To operate the simulator blood control once the system has been powered on you will use the supplied key fob. For upper body bleeding depress once the upper body bleeding button which is the button directly under the fob light. Press again to disengage. The instruction for lower body bleeding is the same but you will use the lower button on the key fob.

4. Replacement of the key fob battery may be necessary over time. This can be done by simply removing the two (2) screws on the back of the key fob and carefully separating the key fob housing. Inside you will find two (2) small lithium powered watch batteries. When replacing batteries be sure to replace with identically numbered batteries.



! " #

Shutting Down the Simulator

Turn off the simulator by pushing in the latching illuminated button on the back of simulator's neck. As soon as you push the button the blue illuminated light will go dark and power will be discontinued.

Removing the Blood Mixture

Blood can be left in the body for short periods of time (up to one week) but should be removed if a longer period of storage is expected.

To remove the blood prior to long-term storage:

1. Turn the simulator on by pushing the latching switch on the neck of the simulator. You will see a blue light when power is applied.
2. Using the key fob controller, turn the bleeding 3' to activate the upper and lower body blood flow until the bleeding stops.

Flushing the System

Should you wish to flush the system:

1. In the hand pump, mix one (1) capful of an antibacterial additive such as Germaben® II-E for every gallon of water.
2. Run this mixture through the simulator until the tank is empty and the water is no longer discharged from the simulator.

Storing after Use

When not in active use, prevent accidental damage by keeping the simulator in the Storage Case. (See *Storing the Simulator* below.)

Performing Interventions

The simulator is designed for training for multiple types of injuries.

Needle Decompression

Supports needle decompression techniques using a 14-gauge needle.



Figure 5: Needle Decompression

When the skin becomes dimpled from multiple punctures, the skin can be repaired. See *Repairing the Skin* below.

Oropharyngeal Airway

Supports oropharyngeal airway techniques using a bag valve mask (BVM) for assisted ventilation.



Figure 6: Oropharyngeal Airway

Nasopharyngeal Airway

A nasopharyngeal airway can be applied to the simulator.



Figure 7: Nasopharyngeal Airway

Wound Packing

Gauze and pressure-dressing materials can be applied to the deep bleeding wounds on the simulator to absorb blood and control bleeding.



Figure 8: Wound Packing

Tourniquet Application

With Bleeding on, standard and junctional tourniquets such as the CAT (Combat Application Tourniquet) and SAM Junctional Tourniquet can be applied to control hemorrhaging at the appropriate sites on the simulator.

Maintaining Your Simulator

Emptying/Cleaning the Bleeding System

Periodically flush the bleeding system by activating the blood flow until the bleeding stops and then infusing the body with a mixture of water and Germaben® II-E (one capful for every gallon of water) to prevent algal/fungal growth.

Repairing the Skin

Over time, the sites for needle decompression will become dimpled.

To repair:

1. Stretch the skin either by pulling gently from the side or by inserting a finger into the humeral or sternal I/O access slits and pushing gently away from the body to spread open the holes.
2. Squeeze a small amount of filler mixture onto the holes and rub in a circular motion to fill the holes.
3. Release the skin and wipe off the excess.



Figure 9: Dimpled Skin from Needle Punctures



Figure 10: Stretching the Punctured Skin



Figure 11: Applying Filler



Figure 12: Rubbing Filler into the Holes

Charging/Replacing the Battery Pack

Simulators are powered by a replaceable, rechargeable, silicone-encased battery pack. For temperature management, the simulators use trickle charge batteries with temperature sensors.

The battery pack is located in a concealed compartment in the inner, upper left thigh. For access, find the seam and slide the skin cover up or down to reveal the compartment. Remove the cover to access the battery

Removing and Installing the Battery Pack

The battery pack should be removed from the simulator during storage and recharging.

To remove the battery pack move the skin cover on the upper left thigh to open and pull the pack from the compartment to access the connector. Detach the connector and remove the pack.

To install the battery pack, reverse this process, moving the skin cover on the upper left thigh to access and open the compartment. Then, attach the pack to the connector, and after resetting the pack in the compartment, close the compartment and replace the skin.

Charging the Battery Pack

Battery packs should only be charged with the charger included with the simulator or otherwise provided by OEI. It takes 3.5 hours to fully charge the trickle charge batteries. The batteries should not remain connected to the charger longer than 14 hours.

To charge a battery pack:

1. Place the charger on a flat surface.
2. Remove the battery pack from the simulator.
3. Connect the battery pack to the charger.
4. Plug into a standard electrical outlet. The red light indicates charging.
5. Charge the batteries for 3.5 to 4 hours (or overnight). The indicator light turns green when the batteries are fully charged and the charger goes into a dormant mode to keep from overcharging.
6. Disconnect the batteries from the charger when fully charged.

Storing the Simulator

Simulators are delivered in custom, wheeled, stackable storage cases.

Storage and Transportation Instructions

- Simulators should be dry, fluid free, and without clothes during long-term storage.
- Battery packs should be removed from the simulator. See *Removing and Installing the Battery Pack* above.
- The straps should be applied by NOT TIGHTENED as this can wear down the moulaged areas.
- Secure all latches.
- The storage cases should NOT BE STACKED more than TWO high.
- Simulators should be stored in a dry environment.

Using the Storage Base for Transport

Cases are inverted with a gurney-like base to field-carry or classroom-carry the simulator. Simulators must be secured with the Velcro® straps during transport.



Figure 13: Gurney-Like Storage Base Used for Transport