

POLICE LINE DO NOT CROSS

Tactical Hemorrhage Control Trainer

The Tactical Hemorrhage Control Trainer is a full-sized, remotely-activated simulator designed for point-of-injury, tactical medicine training for law enforcement and first responders. Ultra-realistic and fully mobile, it is ideally suited for training rapid assessment and treatment of trauma injuries in disaster drill and active shooter simulation scenarios.



Full-body Mannequin

- > Anatomically-accurate soft tissue
- > Ruggedized and water-resistant
- > Remotely-activated pulsatile bleeding
- > Multiple gunshot wounds



Rapid Assessment

- > Airway and breathing
- > Oral and nasal
- > Bleeding
- > Pulses

Airway Management

- > Oropharyngeal and Nasopharyngeal airway



Hemorrhage Control

- > Gunshot wound to the chest with exit wound
 - Chest Seal
 - Needle decompression (mid-clav/mid-ax)
- > Junctional and standard tourniquet application
- > Wound packing